

Release Notes include New Features, Enhancements, Fixes, and Known Issues.

Internal users can use the ISSUE code to view the original Issue in JIRA.

External users can use the menu path in the DESCRIPTION to navigate to the page of the software application where the change was applied. Consult with your Solution Consultant for more information.

NEW FEATURES	
ISSUE	DESCRIPTION OF THE NEW FEATURE
MM-1994	<p>Treatment Plan Editor. Use existing library of treatment plans, modify those or create your own treatment plans from scratch, whatever works best for your customized care program. Connects diagnosis with symptoms, goals, objectives and interventions – all of which are customizable. Systematic and consistent view / edit / print / sign all treatment plans, both on desktop and on iPad.</p> <p>Benefit: Ability for individual treatment options, starting with all substance abuse treatment plans from the DSM-V handbook.</p>
MM-2225	<p>Custom Forms Modernization. Forms can adapt to customer workflow; used throughout the system (InTake, VOB, ASAM, Case Management etc.). Consistent view / edit / sign / print. Same forms can be used on desktop and iPad.</p> <p>Benefit: Ability to tailor forms to specific workflows. Systematic documentation of all aspects of client care is important for care outcomes and compliance.</p>
MM-2441	<p>UI Refresh. Revamp of our UI across three areas (iPad, browser, website), including new logo and color scheme.</p> <p>Benefit: Ease of use through clean, modern and consistent look across all touchpoints.</p>

ENHANCEMENTS	
ISSUE	DESCRIPTION OF THE ENHANCEMENT
MM-0025	Print does not display all items in form
MM-0063	Printing - Logos on group notes
MM-0480	Signature Pad - add support for custom forms
MM-0629	Forms closing upon signature
MM-0691	Window Collapse/Expand Feature
MM-1902	Forms-Printing layout closer to view and printing check boxes
MM-2443	<p>iPad Update.</p> <ul style="list-style-type: none"> - More functionality including all custom forms and treatment plans. - Leverage of iPad specific capabilities (client photos, signature with touch, i.e. no signature pad required, speech recognition). - More consistency with browser-based version (color-coded client cards, filters). - More modern UI theme and alignment with new logo / color scheme. <p>Benefit: More workflows supported on mobile device, saves time and money.</p>

FIXES	
ISSUE	DESCRIPTION OF THE FIX
MM-0733	blank date fields display as 12/31/1969
MM-0797	forms - required "provider dropdown" field - error message
MM-1036	Forms-Tabbing out of a section brings all info back
MM-1575	Client Card - Forms - Error popup is not showing required fields
MM-1576	Client Card - Forms - Form is completed even if required fields are not added
MM-1585	Client Card - Form - Next Back should not be visible for single tabbed form
MM-1605	Client Card - Assessment Form - Missing fields in PRINT
MM-1830	Diagnosis favorites don't show up in the Diagnosis Data List in a custom form
MM-1834	Custom Group Form - Default 'TRUE' property value is not getting saved
MM-2178	By default 'True' control should be checked in encounter view
MM-2179	Forms Radio button required property issue in patient chart
MM-2221	MM-611 Copy Section Button allows multiple unlabeled popup-dialog forms
MM-2222	MM-611 CTRL-Click in Copy Section dialog causes error
MM-2270	Update copyright and improve About dialog
MM-2463	Labs/Imaging Tab - Data not saved in form
MM-2559	Export Patient Export failed and generates error for treatment plan
MM-2572	MAR: Filter menus don't display all items after you select a different entry
MM-2591	VOB Form Comment date is not set with current date

KNOWN ISSUES	
ISSUE	DESCRIPTION OF THE KNOWN ISSUE
MM-0496	Update Policy Holder Information on Insurance Tab
MM-0822	Groups - delete client - group note still appears in chart
MM-0884	combined schedule - Users are still able to overbook 2 groups on 1 schedule
MM-0982	Schedule - group session- Attendees not showing up in recurrence
MM-1084	schedule - dragging an appointment over a conflicting time appointment
MM-1970	VOB form - New Task dialog should not Close without Asking to Save
MM-2104	Appointments - admitting a patient creates a hidden 30 minute appointment for the user

KNOWN ISSUES	
ISSUE	DESCRIPTION OF THE KNOWN ISSUE
MM-2340	Custom Form Builder- Not prompting to save work
MM-2510	iPad: View Scan button in Past Medical History dialog not functional
MM-2513	iPad: Referrer dropdown on Client Info when creating New Patient is hidden by keyboard
MM-2514	iPad: Edit Alert Comment text box content does not move up when keyboard is displayed
MM-2516	iPad: Tx Plan does not save after all Tabs are filled due to unhandled exception
MM-2518	MM-1994 Treatment Diagnosis Code Service exception on selecting favorite diagnosis
MM-2520	MM-1994 Objective data swaps in Overview tab
MM-2521	MM-1994 Validation bypassed on clicking TAB directly
MM-2523	MM-1994 New custom objectives/symptoms/goals/interventions in TPx placement issue
MM-2526	MM-1994 Form is completed even without Signature
MM-2530	MM-1994 Admin Treatment Diagnosis can't be re-select if removed before
MM-2531	MM-1994 Objective Skipping objective selection shows double toast
MM-2539	MM-1994 Admin Exception is thrown in Treatment Diagnosis and Category window
MM-2544	MM-1994 Master TPx TPx form instance status not updated instantly nor deleted in MTx window
MM-2549	iPad: Demographics tab does not have associated Actions menu
MM-2558	Tx Plan: Treatment Plan associated with Diagnosis does not have a name
MM-2562	Safari Signature Date is not fetched correctly
MM-2580	New Custom Form Add Field dialog is opened too low and cut off at bottom
MM-2600	MM-1994 iPad Prepared by is blank in Form list
MM-2603	Old iPad app has issues with new 3.0 server version
MM-2615	iPad not showing Encounters
MM-2623	Forgot Password throwing exception after providing secret answer
MM-2624	RESET password links are not expired
MM-2625	2FA Change password cleared remember device regardless of user check